



CONTROL SYSTEMS

GAS CONTROL AND CLOSING BOX

For main medical gas pipeline connection with the operating rooms, intensive care and patient wards



FEATURES:

- Connection of main pipeline with the consumer
- Consumer lines cut off from the main pipelines
- Gas pressure control
- Sound and light alarm system
- Transmission of alarm signal to central control station
- Lock with emergency opening
- Two models: surface and hidden installation

M-KP24

Technical data:

Characteristics	Model				
	M-KP21	M-KP22	M-KP23	M-KP24	M-KP25
Number of pipelines	1	2	3	4	5
Maximum gas flow in the line, l/min	500	500	500	500	500
Pressure control range, bar	0 ÷ 10	0 ÷ 10	0 ÷ 10	0 ÷ 10	0 ÷ 10
Vacuum control range, bar	-1 ÷ 0	-1 ÷ 0	-1 ÷ 0	-1 ÷ 0	-1 ÷ 0
Alarm signals:					
Light	+	+	+	+	+
Sound	+	+	+	+	+
Medical gas output diameter, mm	15	15	15	15	15
Vacuum output diameter, mm	22	22	22	22	22
Dimensions, L x H x W mm	300x300x114	300x300x114	434x300x114	434x300x114	530x300x114
Weight, kg			9,7		12,8

FEATURES

- Removable frame with doors for easier installation
- Convenient connection of pipelines
- Modification for hidden installation pipelines

M-KP34



- The lock with emergency access shut off valves
- Hidden installation
- Electric section is isolated from gas section

MODELS:

Surface installation

M-KP21
M-KP22
M-KP23
M-KP24
M-KP25
M-KP26


Hidden installation

M-KP31
M-KP32
M-KP33
M-KP34
M-KP35
M-KP36

O₂, AIR5, N₂O, CO₂, VAC, AIR8

- Can be various medical gas pipeline combinations

Convenient installation



2 YEARS GUARANTEE for all gas control and closing boxes

Email: mtlbi@medtech.lt