

## The Most Cost-Effective Solution For Treating Your Medical Waste.



**Decontamination: Proven 8 Log<sub>10</sub> reduction of Bacillus Atrophaeus Spores**

The Envirotek250 produces NO odor. Uses NO steam. NO heat. NO pressure. NO microwaves. It requires NO HEPA filter. NO water softener. NO technical labor to operate.

• **Innovative New Product of 2016** •




## Specifications



Efficacy.....	Exceed STAATT Level III
Cycle Time.....	10 Minutes
Weight.....	1000 kg
Installation.....	Supplied as a mobile unit on wheels
Water Consumption.....	30 Liter per cycle
Electrical Consumption.....	0.35 kwh per cycle
Chamber Capacity.....	250 Liter
Electric Service.....	480 VAC 3 phase–North America 400 VAC 3 phase–Europe
Inlet Water/Drain.....	2 cm Domestic cold water supply.



Envirotek250  
Medical Waste Treatment System  
**Item # ET250**



Envirotek  
Bin Tipper  
**Item # BT250**



E-Cide  
Disinfectant Tablets  
**Item # EC150**



**Corporate Office:**  
**U.S. and Latin America**  
2701-A Hutchison McDonald Road,  
Charlotte, North Carolina 28269  
[info@envirotekwww.com](mailto:info@envirotekwww.com)

**Middle East & North Africa:**  
Office No. 208, RAK Building  
Near AlTwar Center, Al Nahda Rd., Al Qusais 2  
Dubai, UAE  
[info.mena@envirotekwww.com](mailto:info.mena@envirotekwww.com)

**Europe:**  
Cervellati, 3C 40122  
Bologna, Italy  
[info.europe@envirotekwww.com](mailto:info.europe@envirotekwww.com)