



FLOWLIFE™ SUPPLY SOURCES

MEDICAL AIR PLANTS

For centralized compressed air generation in hospitals: for patient's breathing, for surgical instruments, for other medical equipment.



SET OPTIONS:

- 2 or 3 screw type long life compressors
- Pressure receivers are galvanised from inside and outside
- Refrigerating air dryers
- Filtration and pressure reducing unit
- Control unit with integrated alarm system

Technical data:

Model	Voltage, V	Productivity, ltr./min.	Receiver, ltr.	Dew point, °C	Filtration	Pressure to a consumer, bar
M-KS01	230/400V, 50Hz, ~3	2 (3) x 320	500	+3	0,1 µm filter + 0,01 µm filter + activated charcoal	5, 7-8
M-KS02	230/400V, 50Hz, ~3	2 (3) x 460	500	+3		
M-KS03	230/400V, 50Hz, ~3	2 (3) x 650	750	+3		
M-KS04	230/400V, 50Hz, ~3	2 (3) x 950	750	+3		
M-KS05	230/400V, 50Hz, ~3	2 (3) x 1500	900	+3		
M-KS06	230/400V, 50Hz, ~3	2 (3) x 1960	900	+3		
M-KS07	230/400V, 50Hz, ~3	2 (3) x 2620	900	+3		
M-KS08	230/400V, 50Hz, ~3	2 (3) x 3120	2x900	+3		
M-KS09	230/400V, 50Hz, ~3	2 (3) x 4120	2x900	+3		

FILTRATION AND PRESSURE REDUCTION UNIT FOR COMPRESSED AIR

- The set of filters with differential pressure manometers:
 - 0,1 μm
 - 0,01 μm
 - charcoal
 - sterile
- Filter model depending on station's productivity
- Pressure reducers:
 - 500 kPa
 - 800 kPa
- Shut off valves
- Connection to the pipelines diameter $\varnothing 15 - \varnothing 32$



MEDICAL AIR PLANT CONTROL UNIT



- Programmed for control of two or three pumps
- Digital pressure manometer
- Alarm system
- Possibility of connection to central control system
- Dimensions: 444 x 344 x 150 mm

2 YEARS GUARANTEE for all compressed air supply units

Email: mtlbi@medtech.lt